## **SLM-2ER Slit Lamp**

Powerful SLM-2ER professional slit lamp microscope

SLM-2ER professional slit lamp microscope, with a large field of view, good scalability characteristics

- OLarge view range, high magnification, good scalability
- ©Five-step drum zoom design, easier to use
- OMaximum 14mm spot checks greater extent of disease
- ©21.5mm high eye point more comfortable to check
- All optical lens are made from import materials, carried out moisture proof, mould proof and coating process, which prove its good optical performance in different weather condition
- OHigh quality micro-positioning platform, easy to operate
- ©Reserved digital upgrade interface, you can upgraded to a digital slit lamp



Optional teaching tub



Model	SLM-2ER
Microscope Type	Galileo Parallel
Magnification Change Way	Drum Five Magnifications
Eyepiece Magnification	12.5X
Fotal Magnifications	6X,10X,16X,25X,40X
Diopter Adjustment	−5D ~+5D
llumination	Three Steps ( Stepless Optional)
Slit Width	(0-14MM Continuous)
Slit Height	(1-14MM Continuous)
Slit Angle	0° − 180° Adjustable
Slit Inclination Angle	5° 10° 15° 20°
_ight Source	
ight Spot Diameter	0.2mm,1mm,2mm,5mm,10mm,14mm
Filter	Heat Absorption; Grey; Redfree; Cobalt Blue
Fixation	Red LED
Electrical	
llumination Bulb	German OSRAM Halogen Tungsten Lamp LED lamp(optional)
Input Voltage	110V/220V (±10%)
Input Power	60VA
Diffuse Lens	Yes
Preset Lens	Yes
aser Interface	Optional
Applanation Tonometer Interface	Yes
Base movement:	
<b>&lt;</b> :	100mm±3
<i>(</i> :	75mm±3
<u>7</u> :	30mn±3